

## 2019-nCoV Spike S2 ECD Protein

---



<b>Model</b>	STJ11101220
<b>Host</b>	2019-nCoV
<b>Applications</b>	ELISA
<b>Gene ID</b>	43740568
<b>Gene Symbol</b>	S2 ECD
<b>Dilution range</b>	Lyophilized from a 0.22 µm filtered solution of sterile 20 mM PB, 300 mM NaCl, pH 7.0
<b>Note</b>	FOR RESEARCH USE ONLY (RUO).
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store the lyophilized protein at -20°C to -80°C for long term. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.
<b>Alternative Names</b>	S2 protein, Spike glycoprotein Subunit2, S glycoprotein Subunit2, Spike protein S2, COVID-19