

2019-nCoV 3C-like Proteinase (E. coli, N-6His)



Model	STJ61000
Host	E. coli
Applications	ELISA
Immunogen Region	N-6His
Specificity	Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Note	FOR RESEARCH USE ONLY (RUO)
Formulation	Supplied as a 0.2 ml filtered solution of PBS, 1mM EDTA, 10% Glycerol, pH7.4.
Storage Instruction	Store at -20°C. Stable for 6 months. Minimize freeze-thaw cycles.
Alternative Names	M Proteinase3CL Proteinase