
2019-nCoV S Protein RBD (Mammalian, C-6His)

**Description**

2019-nCoV S protein RBD Protein is produced in a mammalian expression system and the target gene encoding Asn331-Val524 is expressed with a 6His tag at the C-terminus.

Model STJ61010

Host Human

Applications ELISA

Immunogen Region Asn331-Val524 , C-6His

Purification Affinity purification chromatography. Greater than 95% as determined by reducing SDS-PAGE.

Note FOR RESEARCH USE ONLY (RUO)

Molecular Weight Mol Mass: 22. 6 kDa, Apparent Mol Mass: 28-35 kDa, reducing conditions.

Formulation Supplied as a 0. 2M filtered solution of PBS, pH7.4.

Storage Instruction Reconstituted protein solution should be stored at -20°C.

Alternative Names 2019-nCov RBD Protein2019-nCoV Spike RBD Protein