

2019-nCoV S Protein NTD (Mammalian, C-6His)



Description	2019-nCoV S Protein NTD is produced in a mammalian expression system with a 6His tag at the C-terminus.
Model	STJ61019
Host	Human
Applications	ELISA
Immunogen Region	C-6His
Purification	Affinity purification chromatography. Greater than 95% as determined by reducing SDS-PAGE.
Note	FOR RESEARCH USE ONLY (RUO)
Molecular Weight	Mol Mass: 32. 4 kDa, Apparent Mol Mass: 40-70 kDa, reducing conditions.
Formulation	Supplied as a 0. 2M filtered solution of PBS, pH7.4.
Storage Instruction	Reconstituted protein solution should be stored at -20°C.
Alternative Names	NA